

How to use:

1. Open up the bags and take out the LED light and mouth tray.
2. Mount the mouth tray onto the front of the lamp.
3. Open the battery lid on the bottom and remove the plastic film in order for the battery to become activated. Both batteries need to be removed to access the plastic.
4. To start LED light press button. Light will automatically go out after 15 minutes.

Replacement Batteries:

2 x CR2032

Teeth Whitening LED Booster Light

Beconfident® - Oral care innovation

Group Beconfident Europe AB
Tennisvägen 31
SE-269 33 BÅSTAD, Sweden

beconfident.com
support@beconfident.com

b® and **beconfident**® are registered trademarks of Group Beconfident Europe AB Designed & developed in EU

